



Ultra-high performance heat insulation transportation box

BioBox PLUS



**Transportation Package
System for
Temperature Sensitive
Medical Substances**



SUGIYAMA-GEN Co.,Ltd.



Sugiyama-Gen Co.,Ltd.

- **Establishment** : 1932
- **Scope of business** :
Medical devices and ancillaries for
healthcare market
- **Main product line**
 - Micro lab. related
 - Animal lab. related
 - Environmental analysis related
 - Glassware for physical/chemical analysis
 - Stainless steel container



Recent Business Focus

- Packing solution for safe and temperature controlled transportation of pharmaceuticals and specimen
 - **BioBox Plus (patent applied)**: transportation system for temperature sensitive medical substances.
examples are as follows;
 - Vaccines
 - General temperature/Freezing sensitive pharmaceuticals
 - Pharmaceuticals/Specimen for clinical trial
 - Specimen for microbiological analysis
 - Others which require strict temperature control
 - **BARRIAPOUCH (with patent)**: pouch for safe transportation of bio-hazard substances



BioBox Plus Presentation

- **Contents:**
 - Requirements about transportation of temperature sensitive medical substances
 - Product configuration
 - Illustration of BioBox Plus
 - Product specifications and features
 - Heat insulation technology
 - Illustration of outer box panel
 - Temperature control technology
 - Temperature control System & Procedure
 - Case study
 - Summary
 - Validation reports on temperature control



Requirements for Transportation of Temperature Sensitive Medical Substances

- **Safe and effective transportation**
- **Stable and homogeneous temperature control, temperature ranging from**
 - $2^{\circ}\text{C} - 8^{\circ}\text{C}$
 - $15^{\circ}\text{C} - 25^{\circ}\text{C}$

(This presentation covers mostly 2-8°C)
- **Durable packaging/storage for long term use**
- **Environmentally friendly system, materials**



Product Configuration

High performance insulation outer box to contain

- (heat sensitive substance)**
- phase change materials (PCM)**
 - cooling agent (coolant/PCM)
 - heating agent (PCM)

-inner box

with pocket for invoice, documents, etc. on surface of lid and front of box

Options:

- Temperature monitoring instrument**
 - Battery driven heater**
 - Anti-freezing agent**
- Note: Outer box itself is widely used in Japan due to superlative heat insulation.



Illustration of **BioBox PLUS** Appearance





BioBox PLUS Components



Inner box
PLUS Inner Box



PCM A : Phase change material
Coolant



PCM B : Phase change material
ThermoPack PLUS 4



Product Specifications

- Name: **BioBox Plus**
- Made of Aluminum with high thermal conductivity
- Patent applied
- Product code: **SBE-P45**
- Dimensions: **STD size**
(customization is available based on customer requirements)
- Outside box**
 - External (mm): 512(+/-10)*388(+/-10)*454(+/-10)height
 - Internal (mm): 418(+/-10)*294(+/-10)*364(+/-10)height
 - (Storage volume (without Inside box): 44.7L
- Inside box**
 - Internal (mm): 372(+/-2)*272(+/-2)*184(+/-2)height
- **Storage volume** (with inside box and PCM): **18.6L**
- **Net weight (Kg): about 5.5Kg plus weight of PCM**
(depending on selection)
- **Materials:** refer to "**Illustration of Outer Box Panel Cross Section**"



Product Features

Outer box: (Heat insulation)

- Five faces and lid are made of superlative vacuum insulation panel and Styrene form layers.
- All corners are covered by the same insulation.
- Lid fits into box with captive labyrinth for air tightness and strength
- Relatively thin thickness of panels

Phase change materials(PCM) : (Temperature control)

- 2 kinds of PCM are employed for cooling and heating
(In case of temperature ranges of 2 °C -8°C)
(Other kind of PCM is for temperature ranges of 15 °C - 25°C)

Inner box: (Heat distribution)

- Inner Box is made of Aluminum plates which bring homogeneity of temperature distribution within inner box




Heat Insulation Technology

Outer Box Panel

- **Heat insulation panel of the outer box**
 - Panel is made of 3 layers.
 - Outer layer: Extruded Polystyrene form
 - Middle layer: Vacuum insulation panel
 - Inner layer: Extruded Polystyrene form



Illustration of Outer Box Panel Cross Section





Temperature Control Technology

- Between the inside of the outer box and the outside of the inner box, 2 kinds of cooling and heating agents(PCM) are put to minimize effect of ambient temperature and effect of cooling agent from heat up and freezing respectively. Both utilize latent heat.
- Appropriate PCMs and their volumes must be selected and associated pre-conditioning is necessary to utilize the phase change from solid to liquid and vice versa.
- Inner box is made of the Aluminum plates.
It can bring even temperature distribution due to good heat energy transmission of Aluminum.



Illustration of BioBox Plus and Components & Procedure

Controlled temperature range (2-8°C)


<Preconditioning>

PCM A : Stabilize PCM in refrigerator at temperature not higher than -15°C over 10 hours

PCM B : Stabilize PCM in a refrigerator at temperature not lower than 4°C over 6 hours

<Procedure>

- ① Install PCM A
- ② Install PCM B
- ③ Put inner box and store heat sensitive substances
- ④ Install PCM B
- ⑤ Install PCM A





One Case Study

- **Assumptions**

- Ambient temperature: 35°C
- Targeted control temperature: 2 °C - 8°C

- **Heat energy transmission**

- Heat energy is transferred from the ambient air to the inside of the box through the insulation panel.
- Incoming heat energy of PCM A becomes latent energy resulting in melting at 0°C.
- PCM A takes heat energy of PCM B resulting in solidifying of PCM B at 3.5°C.




Application Matrix (2-8°C 48Hr)

PCM & Volume VS Ambient Temp.

Control Target Temp. 2-8°C Storage Span 48Hr													
			Ambient temp.(°C)										
Heat control agent			-10	-5	0	5	10	15	20	25	30	35	
Name	Position	Weight (Kg)											
Thermopack Plus4	Top	2											
	Bottom	2											
Cooling agent	Top	2											
	Bottom	2											
Thermopack Plus4	Top	2											
	Bottom	2											
Cooling agent	Top	0.5											
	Bottom	0.5											
Thermopack Plus4	Top	4											
	Bottom	4											
Cooling agent	Top	1											
	Bottom	1											




2-8°C 48 Hr Ambient -10°C






2-8°C **72Hr** **Ambient** **35°C**






2-6°C 240Hr Ambient 25°C





2-8°C 48 Hr Temp fluctuation (-10°C-- +20°C)

2011.01.29~01.31
Testing Agency: SUGIYAMA-GEN Co., Ltd.




BioBox PLUS 2-8°C 72Hr Ambient


2012.12.21~12.25
Testing Agency : SUGIYAMA-GEN Co.,

<Test Condition>

Cold Box : BioBox PLUS
 coolant : Grade 0 °C 500 g × 8 sheets ...①
 inject condition : freezer → in a moment inject BioBox PLUS
 phase change material (Thermo Pack 4) : Grade 4 °C 1500 g × 1 sheet ...②
 inject condition : keeping refrigerator (over 5 °C) → in a moment inject BioBox PLUS
 phase change material (Thermo Pack 4) : Grade 4 °C 2000 g × 1 sheet ...③
 inject condition : keeping refrigerator (over 5 °C) → in a moment inject BioBox PLUS
 coolant : Grade 0 °C 500 g × 12 sheets ...④
 inject condition : freezer → in a moment inject BioBox PLUS



- in the central part of the inner Box
- center of the inside of the bottom of the inner Box
- middle of medial surface of the inner Box
- lower of medial surface of the inner Box
- Sample
- ambient temperature




APEX INTERNATIONAL AICB-444L 2-8°C 96Hr Ambient

2012.12.21~12.25
Testing Agency: SUGIYAMA-GEN Co.,

<Test Condition>

Cold Box : APEX INTERNATIONAL 444-L
Coolant : Cool Pack × 44 sheets (Desinated position of the maker)
inject Condition : keeping refrigerator (Around 3~5°C) → in a moment inject Cold





Summary of **BioBox PLUS**

- **Superlative heat insulation outer box** brought by about 100% box coverage of vacuum insulation panel with durability. This box itself is widely used in Japan due to superlative heat insulation.
- **Effective use of phase change materials** for good temperature control. Preconditioning requires only a normal refrigerator.
- **Aluminum inner box** brings good, homogeneous temperature distribution.